6254421603
BULK DE
**ENGLISH**
DEUTSCH
NEDERLANDS
SUOMI
POLSKI
FRANÇAIS

**Resin-based Dental Restorative Material**

# ESTELITE BULK FILL *Flow*

REF PE01

## ENGLISH

Read all information, precautions and notes before using.

**■PRODUCT DESCRIPTION AND GENERAL INFORMATION**

- ESTELITE BULK FILL Flow is a low viscosity, light-cured, radiopaque composite resin. This low stress flowable material can be placed 4mm increments.
- ESTELITE BULK FILL Flow incorporates Radical-Amplified Photopolymerization initiator technology (RAP technology). RAP technology facilitates a shortened light curing time and ample working time.
- ESTELITE BULK FILL Flow contains 70% by weight (56% by volume) of spherical silica-zirconia filler and composite filler. The composite filler offers reduced polymerization shrinkage stress due to its round shape. The spherical filler contained in ESTELITE BULK FILL Flow (mean particle size: 200nm, particle size range: 100 to 300nm) facilitates excellent gloss retention and wear resistance. ESTELITE BULK FILL Flow contains Bisphenol A di(2-hydroxy propoxy) dimethacrylate (Bis-GMA), Triethylene glycol dimethacrylate (TEGDMA), Bisphenol A polyethoxy methacrylate (Bis-MPEPP), MequinoI, Dibutyl hydroxyl toluene and UV absorber.
- ESTELITE BULK FILL Flow is provided either in SYRINGE or Pre Loaded Tip (PLT).

**■SHADE**

U (Universal), A1, A2, A3, B1

- Please check with your local distributor for available shades.

**■INDICATIONS**

- Direct anterior and posterior restorations
- Cavity lining
- Blocking out cavity undercuts before fabrication indirect restorations
- Repair of porcelain/composite

**■CONTRAINDICATIONS**

ESTELITE BULK FILL Flow contains methacrylic monomers. DO NOT use ESTELITE BULK FILL Flow for patients allergic to or hypersensitive to methacrylic and related monomers or to any of the other ingredients.

**■PRECAUTIONS**

- DO NOT use ESTELITE BULK FILL Flow for any purposes other than those listed in these instructions. Use ESTELITE BULK FILL Flow only as directed herein.
- ESTELITE BULK FILL Flow is designed for sale and use by licensed dental care professionals only. It is not designed for sale nor use by non-dental care professionals.
- DO NOT use ESTELITE BULK FILL Flow if the safety seals are broken or appear to have been tampered with.
- If ESTELITE BULK FILL Flow causes an allergic reaction or hypersensitivity, discontinue its use immediately.
- Use examination gloves (plastic, vinyl, or latex) at all times when handling ESTELITE BULK FILL Flow to avoid the possibility of allergic reactions from methacrylic monomers. Note: Certain substances/materials may penetrate through examination gloves. If ESTELITE BULK FILL Flow comes in contact with the examination gloves, remove and dispose of the gloves, and wash hands thoroughly with water as soon as possible.
- Avoid contact of ESTELITE BULK FILL Flow with eyes, mucosal membrane, skin and clothing.
  - If ESTELITE BULK FILL Flow comes in contact with the eyes, thoroughly flush eyes with water and immediately contact an ophthalmologist.
  - If ESTELITE BULK FILL Flow comes in contact with the mucosal membrane, wipe the affected area immediately, and thoroughly flush with water after the restoration is completed.
  - If ESTELITE BULK FILL Flow comes into contact with the skin or clothing, immediately saturate the area with an alcohol soaked cotton swab or gauze.
  - Instruct the patient to rinse his mouth immediately after treatment.
- ESTELITE BULK FILL Flow should not be ingested or aspirated. Ingestion or aspiration may cause serious injury.
- To avoid the unintentional ingestion of ESTELITE BULK FILL Flow, do not leave it unsupervised within the reach of patients and children.
- Clean the placement instruments and brushes with alcohol after use.
- When using a light-curing unit, protective eye shields, glasses or goggles should be worn at all times.
- ESTELITE BULK FILL Flow is available in semi-translucent shades. If applied in cavity Class III or IV direct restorations, the cured ESTELITE BULK FILL Flow may appear dark. Avoid using ESTELITE BULK FILL Flow for the direct restoration of Class III or IV when opacity is paramount.
- Light-cure ESTELITE BULK FILL Flow extra orally to confirm the time needed for complete hardening of ESTELITE BULK FILL Flow with your light-curing unit prior to doing the restoration.
- This product is designed to be used at room temperature (18 - 30<sup>°</sup> C / 62 - 84<sup>°</sup> F). Allow product to reach room temperature prior to use. Cold material may be difficult to extrude.

- CURING TIME AND DEPTH**

ESTELITE BULK FILL Flow enables a 4mm depth of cure (U shade) by light-curing for 10 seconds with Halogen or LED light-curing unit (Light intensity: 800mW/cm<sup>2</sup> or more, wavelength: 400-500nm). The depth of cure was determined based on tests performed in accordance with section 7.10 “depth of cure” of ISO 4049:2009.
- SPECIAL NOTES FOR THE USE OF PLT**

1) PLTs are designed for single patient use only. Do not re-cap and/or re-use the PLT tip once material has been dispensed for that patient.

2) Dispensers are not provided for ESTELITE BULK FILL Flow. Use a dispenser that fits the PLT of ESTELITE BULK FILL Flow. For proper use and proper disinfection, see manufacturer’s instructions.

3) Use light, controlled pressure to prevent any continued extrusion of material following pressure release.

**■STORAGE**

- Store ESTELITE BULK FILL Flow at a temperature between 0 - 25°C (32 - 77°F).
- AVOID direct exposure to light and heat.
- DO NOT use ESTELITE BULK FILL Flow after the indicated date of expiration on the syringe or PLT package.

**■DISPOSAL**

To safely dispose of excess ESTELITE BULK FILL Flow, extrude unused portion from SYRINGE or PLT and light-cure before disposal.

**■CLINICAL PROCEDURE**

**1. Cleaning**

Thoroughly clean the tooth surface with a rubber cup and a fluoride-free paste then rinse with water.

**2. Shade Selection**

Prior to isolation of tooth, select the appropriate shade of ESTELITE BULK FILL Flow using shade card provided for ESTELITE BULK FILL Flow.

**3. Isolation**

A rubber dam is the preferred method of isolation.

**4. Cavity Preparation**

Prepare the cavity and rinse with water.

**5. Pulp Protection**

Calcium hydroxide should be applied if the cavity is in close proximity to the pulp. DO NOT use EUGENOL-BASED MATERIALS to protect the pulp as these materials will inhibit curing ESTELITE BULK FILL Flow. ESTELITE BULK FILL

Flow is not indicated for direct pulp capping.

**6. Bonding system**

Apply a light-cured bonding system according to its manufacturer’s instructions.

- Do not use self-cured bonding systems. We advise that light-cured bonding systems should be used with ESTELITE BULK FILL Flow since other types of bonding systems such as self-cured or dual-cured bonding systems may not provide effective adhesion when used with light-cured composite resins. If you want to use self-cured or dual-cured bonding systems, please confirm that such bonding systems are compatible with the light-cured composite resins being used.
- In case of porcelain/composite repairs be sure to condition the surface of porcelain/composite with a silane coupling agent BEFORE performing the bonding procedure.

**7. Dispensing**


7-1. PLT

Please read **■SPECIAL NOTES FOR THE USE OF PLT** before use.

- Load the PLT into a dispenser that fits the PLT.
- Remove the PLT cap.

7-2. SYRINGE

- Hold the nozzle and turn the cap counterclockwise for removal. [1]
- Hold the nozzle and place the Dispensing Tip onto it. Turn the Dispensing Tip clockwise until it is securely locked. [2]

- Accessories and Dispensing Tip are exclusively designed for ESTELITE BULK FILL Flow. Avoid using Dispensing Tips from other brands with ESTELITE BULK FILL Flow.
- Dispensing Tips are also sold separately.
- Use gloves when attaching the Dispensing Tip.

**8. Filling and contouring**

Cavity can be directly filled from the Dispensing Tip/PLT. Increments should not exceed 4mm at a time.

- DO NOT apply excessive force to extrude the material from the syringe or PLT immediately after removal from the refrigerator.
- DO NOT mix ESTELITE BULK FILL Flow with other brands of resin composite or other shades of ESTELITE BULK FILL Flow to avoid incomplete cure or entrapment of air bubbles.
- In case of syringe, after extruding the paste, remove and discard the Dispensing Tip, wipe the nozzle with gauze and replace the cap immediately.

**9. Curing**

Light-cure each increment for at least the indicated time (please see the curing recommendations below). Light-cure each layer (up to 4mm increment) using Halogen light-curing or LED light-curing unit with a wavelength of 400-500 nm. Be sure the light-curing unit with an intensity of 800mW/cm<sup>2</sup> or more. Hold the light guide tip as close to the restorative as possible during light exposure.

Shade	Increment Depth	Light Output	Curing Time
All shades	4mm	≥800mW/cm <sup>2</sup>	10 sec.

- If other brands of composite resins are layered over the cured composite, follow their instructions.

**10. Finishing**

Shape and polish the restoration. For finishing, use finishing discs and/or fine finishing diamond points. Use metal finishing strips or vinyl polishing strips for proximal surfaces. For polishing, polish with rubber points or any suitable polishing tools. For final polishing, use felt discs or cotton wheels with polishing paste, or suitable polishing tools.

IMPORTANT NOTE: The manufacturer is not responsible for damage or injury caused by improper use of this product. It is the personal responsibility of the dental professional to ensure the product is suitable for application before use. Specifications are subject to change without notice. When the product specification changes, the instructions and precautions may change also.

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